

Lec 130 Litre Basic Laboratory Refrigerator (LR111)

Product Code: LEC444441207

Brand: Lec Medical

The Lec 130 Litre Basic Laboratory Refrigerator (LR111) has the power to easily slide under a general laboratory bench. Supplied with a door lock and two keys, this freezer provides laboratories with great security allowing only approved users to open it. Finished in white and spark-free, it has adjustable temperature control and can auto defrost. It also contains three shelves and useful drawers.



Features:

- **Spark-Free Interior:** Offers safe and reliable storage of volatile liquids, reagents, chemicals, poisons, and samples with a spark-free interior that reduces the risk of internal explosions.
- **Lock:** Security is essential in laboratories, so this freezer comes complete with a door lock and two keys, ensuring that only authorised individuals are able to access it.
- Adjustable thermostat.
- Auto defrost.

Specifications:

Capacity (litres)	130
Temperature Range (°C)	2 to 8
Class	Climate Class N
Dimensions (w x d x h)	545 x 554 x 845mm
Electrical Supply	230V, 50/60Hz
Colour/Finish	White
Reversible Door?	Yes
Door Lock?	Yes
Drawer Number	1 Plastic
Shelf Number	3 Glass
Handle Type	Recessed
Defrost Capability	Auto Defrost
Annual Energy Consumption (kwh)	117
Thermostat	Adjustable
Spark-Free?	Yes
HFC Free?	Yes
CFC Free?	Yes