

Medline Analogue Hotplates and Stirrers

Brand: Medline Scientific

Medline's Analogue Hotplates and Stirrers offer a large heating area of 180 x 180mm with safe controlled heating up to 400°C.

Constructed with chemical resistant ceramic coated steel, this range of lab hotplates and stirrers deliver excellent accuracy, reproducible results, and have a secure design.



Specifications:

Product Code(s)	MS-100	MS-200	MS-300
Description	Hotplate	Stirrer	Hotplate Stirrer
Case Construction	Ceramic coated steel (chemical resistant)		
Stage Area (mm)	180 x 180		
Temperature Range (°C)	Amb. to 400	-	Amb. to 400
Speed Range (rpm)	-	100 to 1500	
Retort Rod (12.5mm Ø)	400mm/760mm (optional)		
Dimensions (w x d x h)	410 x 210 x 150mm		
Weight (kg)	4		5