

Medline Rotary Vane Vacuum Pump (CCMD4)

Product Code: MED-141C

Brand: Medline Scientific

Equipped with a direct connect oil seal to prevent air reflux, Medline's Rotary Vane Vacuum Pump (CCMD4) is specially designed for vacuum ovens to achieve a higher level of vacuum pressure.

Ideal for low laboratory applications.

An optional oil mist trap is available to prevent oil in-flow. Recommended for use with Jeio Tech's Vacuum Ovens (OV Series).

Specifications:

Stages	2
Pumping Speed (L/min.)	140
Ultimate Vacuum (m/bar)	5×10^{-3}
Gas Ballast	Yes
Motor Power (W)	370
Max. Water Vapour Pumping Rate (kg/hr.)	400
Rotational Motor Speed (rpm)	2880
Oil Capacity (ml)	450
Recommended Oil (litres)	1.1
Inlet/Outlet Port (DN)	8mm Hose Tail
Oil Mist Eliminator (Exhaust Outlet)	Fitted
Power Compatibility (Single Phase)	110/240V Dual
Dimensions (w x l x h)	178 x 375 x 330mm
Weight (kg)	12

